

Product Datasheet

Aloc-D-Ala-Phe-Lys(Aloc)-PAB-PNP

Catalog No. **BP-00044** · CAS **253863-34-2** · Form **free base** · Physical state **solid** · Version **1.0** · Revision **2026-10-08**


FOR RESEARCH USE ONLY. NOT FOR HUMAN OR VETERINARY USE.

Target ADC Linkers	Research area Cancer
------------------------------	--------------------------------

Molecular formula **C₄₀H₄₆N₆O₁₂**

Molecular weight **802.8262**

STRUCTURE



DESCRIPTION

Aloc-D-Ala-Phe-Lys(Aloc)-PAB-PNP is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs).

Tel: +86 20 31130026 · Fax: +86 20 36822386 · tech@biochemprobe.com · www.biochemprobe.com

BiochemProbe is the catalog brand of Guangzhou Yuyuan Biotech Co., Ltd. · For research use only.

Lot CoA and SDS are separate documents.