

Product Datasheet

SGC6870N

Catalog No. **BP-00147** · CAS **2561471-15-4** · Form **free base** · Physical state **solid** · Color **White to off-white** · Version **1.0** · Revision **2026-10-08**


FOR RESEARCH USE ONLY. NOT FOR HUMAN OR VETERINARY USE.

Target Histone Methyltransferase	Research area Cancer
--	--------------------------------

Molecular formula **C₂₃H₂₁BrN₂O₂S**

Molecular weight **469.3940**

STRUCTURE



DESCRIPTION

SGC6870N is inactive against PRMT6 and can be utilized as a negative control. SGC6870N is an inactive enantiomer of SGC6870.

REFERENCES

[1]. Yudao Shen, et al. A First-in-Class, Highly Selective and Cell-Active Allosteric Inhibitor of Protein Arginine Methyltransferase 6. J Med Chem. 2021, 64, 7.

Tel: +86 20 31130026 · Fax: +86 20 36822386 · tech@biochemprobe.com · www.biochemprobe.com

BiochemProbe is the catalog brand of Guangzhou Yuyuan Biotech Co., Ltd. · For research use only.

Lot CoA and SDS are separate documents.