

Product Datasheet

EED ligand 1

Catalog No. **BP-01279** · CAS **2762275-20-5** · Form **free base** · Physical state **solid** · Color **white to off-white** · Version **1.0** · Revision **2026-09-14**


FOR RESEARCH USE ONLY. NOT FOR HUMAN OR VETERINARY USE.

Target Histone Methyltransferase	Research area Cancer
--	--------------------------------

Molecular formula **C₁₉H₁₉FN₈O**

Molecular weight **394.4056**

STRUCTURE



DESCRIPTION

EED ligand 1 is a diverse, potent, and efficacious inhibitor that target the EED subunit of the polycomb repressive complex 2 (PRC2) methyltransferase.

SOLUBILITY

Soluble in DMSO

REFERENCES

[1]. Bagal SK, et al. Diverse, Potent, and Efficacious Inhibitors That Target the EED Subunit of the Polycomb Repressive Complex 2 Methyltransferase. J Med Chem. 2021 Nov 22.

Tel: +86 20 31130026 · Fax: +86 20 36822386 · tech@biochemprobe.com · www.biochemprobe.com

BiochemProbe is the catalog brand of Guangzhou Yuyuan Biotech Co., Ltd. · For research use only.

Lot CoA and SDS are separate documents.