

Product Datasheet

EED226-COOH

Catalog No. **BP-02022** · CAS **2467965-71-3** · Form **free base** · Physical state **solid** · Color **off-white to light yellow** · Version **1.0** · Revision **2026-09-11**


FOR RESEARCH USE ONLY. NOT FOR HUMAN OR VETERINARY USE.

Target Ligands for Target Protein for PROTAC; Histone Methyltransferase	Research area Cancer
---	--------------------------------

Molecular formula **C₁₇H₁₃N₅O₃**

Molecular weight **335.3168**

STRUCTURE



DESCRIPTION

EED226-COOH is an EED226-derived ligand for target protein EED ligand for PROTAC, binds to a ligand for VHL via linker to form UNC6852 to degrade PRC2.

BIOLOGICAL ACTIVITY

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

SOLUBILITY

Soluble in DMSO : 125 mg/mL

REFERENCES

[1]. Potjewyd F, et al. Degradation of Polycomb Repressive Complex 2 with an EED-Targeted Bivalent Chemical Degradar. Cell Chem Biol. 2020 Jan 16;27(1):47-56.e15.

Tel: +86 20 31130026 · Fax: +86 20 36822386 · tech@biochemprobe.com · www.biochemprobe.com

BiochemProbe is the catalog brand of Guangzhou Yuyuan Biotech Co., Ltd. · For research use only.

Lot CoA and SDS are separate documents.