

Product Datasheet

PH-HG-005-5

Catalog No. **BP-02171** · CAS **2810895-93-1** · Form **free base** · Physical state **solid** · Color **Light yellow to yellow** · Version **1.0**
· Revision **2026-09-04**


FOR RESEARCH USE ONLY. NOT FOR HUMAN OR VETERINARY USE.

Target Drug-Linker Conjugates for ADC	Research area Cancer
---	--------------------------------

Molecular formula **C₄₁H₄₇N₇O₁₀**

Molecular weight **797.8528**

STRUCTURE



DESCRIPTION

PH-HG-005-5 (compound 16c) is a derivative of SN-38 and can be used as Drug-Linker Conjugates for ADC. PH-HG-005-5 can conjugate to targeting peptides for ADCs synthesis .

SOLUBILITY

Soluble in DMSO

REFERENCES

[1]. Chen Jian, et al. EphA2 bicyclic peptide ligand and conjugate thereof. World Intellectual Property Organization, WO2023061482 A1. 2023-04-20.

Tel: +86 20 31130026 · Fax: +86 20 36822386 · tech@biochemprobe.com · www.biochemprobe.com

BiochemProbe is the catalog brand of Guangzhou Yuyuan Biotech Co., Ltd. · For research use only.

Lot CoA and SDS are separate documents.