

Product Datasheet

E3 ligase Ligand 18

Catalog No. **BP-02198** · CAS **2241669-88-3** · Form **free base** · Physical state **solid** · Color **White to off-white** · Version **1.0** · Revision **2026-09-07**


FOR RESEARCH USE ONLY. NOT FOR HUMAN OR VETERINARY USE.

Target Ligands for E3 Ligase	Research area Cancer
--	--------------------------------

Molecular formula **C₁₇H₁₇ClN₆O₄S**

Molecular weight **436.8730**

STRUCTURE



DESCRIPTION

E3 ligase Ligand 18 is a ligand for E3 ubiquitin ligase. E3 ligase Ligand 18 can be connected to the ligand for protein by a linker to form PROTACs. PROTACs are inducers of ubiquitination-mediated degradation of cancer-promoting proteins[1].

SOLUBILITY

soluble in DMSO

REFERENCES

[1]. Heterocyclic compound. WO2018143403A1.

Tel: +86 20 31130026 · Fax: +86 20 36822386 · tech@biochemprobe.com · www.biochemprobe.com

BiochemProbe is the catalog brand of Guangzhou Yuyuan Biotech Co., Ltd. · For research use only.

Lot CoA and SDS are separate documents.