

Product Datasheet

E3 ligase Ligand 43

Catalog No. **BP-02391** · CAS **2821795-71-3** · Form **free base** · Physical state **solid** · Version **1.0** · Revision **2026-09-14**


FOR RESEARCH USE ONLY. NOT FOR HUMAN OR VETERINARY USE.

Target Ligands for E3 Ligase	Research area Cancer
--	--------------------------------

Molecular formula **C₂₂H₂₈N₆O₄S**

Molecular weight **472.5610**

STRUCTURE



DESCRIPTION

E3 ligase Ligand 43 serves as the E3 ubiquitin ligase ligand for ACB13. E3 ligase Ligand 43 can be utilized for the synthesis of PROTACs.

SOLUBILITY

Soluble in DMSO

REFERENCES

[1]. Popow J, et al. Targeting cancer with small-molecule pan-KRAS degraders[J]. Science, 2024, 385(6715): 1338-1347.

Tel: +86 20 31130026 · Fax: +86 20 36822386 · tech@biochemprobe.com · www.biochemprobe.com

BiochemProbe is the catalog brand of Guangzhou Yuyuan Biotech Co., Ltd. · For research use only.

Lot CoA and SDS are separate documents.